

Demo Rad-Control Barcode

Instructions | 2021

Note: This is a simple demo tool created for presentation purposes only.

To get started you need:

1. **8-digit barcode** (see attachment) or any other 8-digit international type 128 barcode for simulating the scanning of the blood product
2. One **sample blood bag**
3. Any suitable **barcode scanner** – we recommend a CCD barcode scanner such as the Datalogic TouchTM 65/90 Pro Contact Reader or any similar device.
4. One **irradiated Rad-Control Barcode label**

Prepare for the presentation:

- Place the 8-digit barcode on the sample blood bag.
This is only for demonstration purposes. In real life the blood bank will have a more complex barcode on a blood product.
- Place the irradiated Rad-Control Barcode label on the sample blood bag.
- Plug the Barcode scanner into your computer.
- Open the zip file RadControl_V8
- Open the file “RadControlDemo” and switch to full screen.

NOTE: The cursor needs to be in the program /program window at all times.

Start the Rad-Control Barcode Demo:

- This is your initial screen

Number	Info	Irradiation result	Date	Operator
129-55	CP2D/AS-3 Double Blood Bag	OK	2021-05-31	MGR
129-53	CP2D/AS-3 Double Blood Bag	OK	2021-05-31	MGR
129-52	CP2D/AS-3 Double Blood Bag	OK	2021-05-31	MGR

- Explain the situation you are going to simulate: a blood bag has just been irradiated and taken out of the irradiator. The results will be automatically stored in your data base to allow for a more efficient and streamlined workflow.
- **STEP 1 - documentation of the blood product:**
Scan the 8-digit barcode.
- Your customer now sees the blood product information being entered into the system on the left side of the screen. The irradiation result is still open and needs to be confirmed next.

Number	Info	Irradiation result	Date	Operator
129-62	CP2D/AS-3 Double Blood Bag	open	2021-05-31	MGR
129-55	CP2D/AS-3 Double Blood Bag	OK	2021-05-31	MGR
129-53	CP2D/AS-3 Double Blood Bag	OK	2021-05-31	MGR
129-52	CP2D/AS-3 Double Blood Bag	OK	2021-05-31	MGR

- **STEP 2 - Confirmation of blood irradiation:**
Scan the barcode from the irradiated indicator label on your sample blood bag. The irradiation result changes from “open” to “OK”.

Product

No.: 129-62

Info: CP2D/AS-3 Double B Bag

Date: 2020-12

Produced by: Haemonetics Corporation

Produced for: Haemonetics Manufacturing

SCAN blood product

SCAN indicator

SCAN blood product ok

SCAN indicator OK

CHANGE STATUS

NEXT

Number	Info	Irradiation result	Date	Operator
129-62	CP2D/AS-3 Double Blood Bag	OK	2021-05-31	MGR
129-55	CP2D/AS-3 Double Blood Bag	OK	2021-05-31	MGR
129-53	CP2D/AS-3 Double Blood Bag	OK	2021-05-31	MGR
129-52	CP2D/AS-3 Double Blood Bag	OK	2021-05-31	MGR

- Click on “Change status” to simulate that the status of the blood product can either be changed manually or automatically.

Product

No.: 129-62

Info: CP2D/AS-3 Double B Bag

Date: 2020-12

Produced by: Haemonetics Corporation

Produced for: Haemonetics Manufacturing

SCAN blood product

SCAN indicator

SCAN blood product ok

SCAN indicator OK

CHANGE STATUS

NEXT

Number	Info	Irradiation result	Date	Operator
129-62	CP2D/AS-3 Double Blood Bag	OK	2021-05-31	MGR
129-55	CP2D/AS-3 Double Blood Bag	OK	2021-05-31	MGR
129-53	CP2D/AS-3 Double Blood Bag	OK	2021-05-31	MGR
129-52	CP2D/AS-3 Double Blood Bag	OK	2021-05-31	MGR

- If needed, repeat the demonstration by clicking on “Next”.

Contact us for a **training session** prior to your first demo presentation!

www.onpointmedical.com